




# neodisher® Alka 300



Alkaline detergent for the automated reprocessing of thermostable and thermolabile instruments



Liquid concentrate

### Fields of application:

- Automated cleaning of sterilisation containers, surgical instruments and anaesthetic equipment
- Manual intensive cleaning of surgical instruments in immersion baths
- Manual intensive cleaning of bedpans by wiping with the fleece wipe dispenser system neodisher wipes RTF
- Automated cleaning of cages in animal facilities

### Performance spectrum:

- Removes blood, protein residues and organic residues
- Suitable for instrument steel, chrome nickel steel, chromium-plated and nickel-plated materials, pure aluminium, aluminium-magnesium alloys and commonly used plastics
- Suitability tests should be carried out before use in the case of special instruments, e.g. optical instruments or anaesthetic and respiration materials, as well as anodised aluminium.

### Special properties:

- With active chlorine for intensive cleaning performance
- With oxidative decomposition action for stubborn organic residues
- Usable with any level of water hardness, but generally used with softened or deionised water

### Application and dosage:

neodisher Alka 300 can be used in washer disinfectors. For intensive cleaning application in immersion baths or via wiping is possible.

The following parameters are recommended for the different application possibilities:

Automated cleaning of sterilisation containers, surgical instruments and anaesthetic equipment	2 - 5 ml/l (0.2 - 0.5 %)*
Automated cleaning of cages in animal facilities	2 - 5 ml/l (0.2 - 0.5 %)*
Manual intensive cleaning of surgical instruments and bedpans	For application parameters for the intensive cleaning of surgical instruments and bedpans please contact your neodisher specialist adviser or our Applications Engineering Department neodisher

\*The dosing amount depends among other things on the degree of soiling.

Use suitable dosing devices.

### General instructions for use:

- The neodisher Alka 300 solution has to be rinsed off completely (preferably with deionised water).
- Do not mix with other products.
- Rinse out dosing system including hose with water (preferably deionised) before changing product.
- For professional use only.
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664, the cleaning



# neodisher<sup>®</sup> Alka 300

recommendations given by the manufacturer of the animal cages as well as the recommendations of the working group for cage processing" (AK KAB) in the current issue of the AK KAB brochure "Cage Processing in Animal Facilities properly done".

- The instructions given by the manufacturer of the washer disinfectant are to be observed.

Dispose only when container is empty and closed.  
For disposal of product residues, refer to the Safety Data Sheet.

MB 3103/3-3

Date of issue 05/2021

## Technical data:

Density	approx. 1.4 g/cm <sup>3</sup> (20 °C)
pH-range	11.1 - 11.9 (2 - 5 ml/l, determined in deionised water, 20 °C)
Viscosity	< 50 mPa s (concentrate, 20 °C)
Titration factor	1.04 (in accordance with neodisher titration method)

## Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:


- < 5 % chlorine-based bleaching agents
- 15 - 30 % phosphates

## CE-mark:

neodisher<sup>®</sup> Alka 300 complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

## Storage information:

Always store at a temperature between -10 °C and 25 °C. Keep away from sunlight.  
Usable for 1 year when stored as recommended.  
For expiry date refer to the stamp mark on the label behind the hourglass symbol .

## Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available at [www.drweigert.com](http://www.drweigert.com) under the category "Service/Downloads".

With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.