neodisher[®] BP

pH-neutral rinse aid with corrosion protection and lubrication
for automated reprocessing of bedframes, surgical tables and trolleys
Liquid concentrate

Main fields of application:

Rinsing and lubrication of bedframes, surgical tables and trolleys in decontamination systems in hospitals and nursing homes.

Characteristics:

neodisher BP is a pH-neutral rinse aid and lubricant based on special oils for corrosion protection. neodisher BP ensures perfect drying of bedframes, surgical tables and trolleys. Moving parts such as brakes, Bowden cables and rollers are kept smooth and at the same time lubricated.

Application and dosage:

For decontamination systems using a circulation or fresh water system:

1 - 4 ml/l is dosed into the rinsing water via suitable dosing devices, depending on the process parameters of the decontamination system and the material of the bedframes, surgical tables and trolleys.

Optimum adjustment of the automatic dosing must be determined on an empirical basis.

Do not mix with other products. Rinse out dosing system including suction hose with water and the rinse aid before changing product.

Reprocessing should comply with all ordinances pursuant to the Medical Device Directive and should be performed with appropriate validated processes.

Please observe the instructions and reprocessing recommendations of the manufacturers of the surgical tables, bedframes and trolleys as well as the recommendations given by the "Working Group for Bedframe and Cart Decontamination Systems" (Arbeitskreis Bettgestell- und Wagendekontaminationsanlagen [AK-BWA]) in the current issue of the AK BWA brochure "Automated Decontamination".

The instructions given by the manufacturer of the decontamination system are to be observed.

For professional use only.

Technical data:

Density (20 °C): 0.9 g/cm³

pH-values according to neodisher BP application solutions at 20 °C in a variety of water qualities:

Concentration of neodisher BP	deionised water	softened from 15 °d to 0 °d
0.0 ml/l	5.1	7.6
1.0 ml/l	5.9	7.6
4.0 ml/l	5.9	7.6

Viscosity (concentrate, 20 °C): < 100 mPa s

Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:
> 30 % non-ionic surfactants, aliphatic hydrocarbons

CE-mark:

neodisher BP conforms to the European directive 93/42/EEC, Annex I concerning medical devices.

Storage information:

Always store at a temperature between 0 °C and 30 °C. At temperatures below 10 °C reversible clouding may occur. Usable for 2 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol

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With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.

neodisher[®] BP

pH-neutral rinse aid with corrosion protection and lubrication
for automated reprocessing of bedframes, surgical tables and trolleys
Liquid concentrate

Hazard and precautionary statements:

For safety information see EC Safety Data Sheets. These are available at www.drweigert.com under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

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Chemische Fabrik Dr. Weigert GmbH & Co. KG