




# neodisher<sup>®</sup> FA Dental



## Alkaline detergent for the automated reprocessing of dental instruments

### Liquid concentrate

#### Fields of application:

- Automated cleaning of dental instruments and dental utensils such as trays, cups for mouth rinsing and kidney dishes

#### Performance spectrum:

- Removes blood, protein, bone residues and ointments
- Provides excellent material protection and a high level protection against corrosion with dental instruments and dental utensils made of stainless steel
- Not suitable for the automated cleaning of dental handpieces. For this purpose we recommend the detergent neodisher MediClean Dental

#### Special properties:

- High material protection
- Free of oxidising agents and surfactants

#### Application and dosage:

neodisher FA Dental is used in washer disinfectors for dental instruments.

The dosage depends among other things on the degree of contamination of the items to be washed and is 4 - 5 ml/l.

Use suitable dosing devices for dosing neodisher FA Dental.

Residues of composite, dental cement and compounds for dental imprints must be removed before they set.

For avoiding water stains the use of deionised water in the final rinse is recommended.

The neodisher FA Dental solution has to be rinsed off completely with water (preferably with deionised water).

#### Notes on application:

- For professional use only
- Do not mix with other products
- Rinse out dosing system including suction hose with water before changing product
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes
- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664.
- The instructions given by the manufacturer of the washer disinfectant are to be observed

#### Technical data:

pH value	11.8 - 12.0 (4 - 5 ml/l, determined in deionised water, 20 °C)
Viscosity	< 50 mPa s (concentrate, 20 °C)
Density	approx. 1.5 g/cm <sup>3</sup> (20 °C)
Titration factor	0.68 (in accordance with the neodisher titration method)

#### Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:  
15 - 30 % phosphates

CE-mark:  


neodisher FA Dental complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.



# neodisher<sup>®</sup> FA Dental

## Storage information:

Sensitive to frost below - 15 °C. Always store at a temperature between -15 °C and 30 °C. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

## Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available e.g. at [www.drweigert.com](http://www.drweigert.com) under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

MB 4124/3-2

Date of issue: 02/2021

---

With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.