




neodisher[®] IP Konz



Care product for instruments for use in washer disinfectors



Liquid concentrate

Fields of application:

- Basic care during automated reprocessing of surgical instruments in the final rinse step

Performance spectrum:

- Especially suitable for use in washer disinfectors for lubricating surgical instruments prior to sterilisation
- No persisting germs in the extremely thin oil film
- Not suitable for instruments that are sterilised with hot air

Special properties:

- Based on highly purified paraffin oil and emulsifiers
- Steam sterilisation is not inhibited
- Residues do not accumulate on the steel surfaces and also do not resinify during the sterilisation process
- When neodisher IP Konz is used as intended, the biocompatibility of the treated instruments is not affected. (Biological evaluation of medical devices according to EN ISO 10993-1). The active substances are in compliance with applicable requirements FDA 21 CFR 176.210 (Emulsifier) and 21 CFR 172.878 (white paraffin) and USP/EUP/DAB (white paraffin).

Application and dosage:

neodisher IP Konz is intended for use in washer disinfectors. The following parameters are recommended for automated use of neodisher IP Konz:

Application in the final rinse step	1 - 2 ml/l, at least 50 °C
-------------------------------------	----------------------------

Suitable dosing devices should be used.

During the subsequent cleaning the neodisher IP Konz film is completely removed; therefore neodisher IP Konz must be re-applied with each disinfection and cleaning process.

neodisher IP Konz is not recommended for use in immersion baths because of the risk of uncontrolled germ contamination.

Alternatively, we recommend the use of neodisher IP Spray after automated reprocessing for direct lubrication of hinged instruments.

Notes on application:

- For professional use only!
- Do not mix with other products.
- Rinse out dosing system including suction hose with water and product before changing product.
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Please observe the medical device manufacturer's recommendations for use according to the requirements of DIN EN ISO 17664.
- Please observe the instructions given by the manufacturer of the washer disinfectant.

Technical data:

Density	approx. 0.9 g/cm ³ (20 °C)
Viscosity	< 50 mPa s (concentrate, 20 °C)

Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:
15 – 30 % non-ionic surfactants



neodisher[®] IP Konz


> 30 % aliphatic hydrocarbons (white paraffin oil DAB)

CE-mark: 

neodisher IP Konz complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

Storage information:

Always store at a temperature between 0 °C and 30 °C. When stored below 10 °C reversible clouding may occur. Usable for 2 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available e.g. at www.drweigert.com under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

MB 4306/03-2 Revision Date 04/2020

Die Angaben dieses Merkblattes basieren auf unseren derzeitigen Kenntnissen und Erfahrungen. Sie befreien den Verwender nicht von eigenen Prüfungen und Versuchen. Eine rechtlich verbindliche Zusicherung bestimmter Eigenschaften kann hieraus nicht abgeleitet werden.