




# neodisher® MediDry special



Rinse aid for the automated reprocessing of thermostable and thermolabile instruments



## Liquid concentrate

### Field of application:

Rinsing of medical devices such as surgical instruments incl. ophthalmological instruments, MIS instruments, anaesthesia utensils, instrument containers, implants as well as baby bottles in washer disinfectors

### Performance spectrum:

- Supports the drying, especially of plastic materials
- Suitable for all materials, especially for accessories of anaesthesia utensils made of polysulfone (PSU) and polyphenylsulfone (PPSU) such as connectors of laryngeal masks, adaptors and container lids
- When neodisher MediDry special is used as intended the biocompatibility of the treated instruments is not affected (biological assessment of medical devices according to EN ISO 10993-1)

### Special properties:

- Free of surfactants
- Acidic, neutralises alkaline residues
- Reduces the formation of stains when softened water is used in the final rinse
- Excellent material compatibility; prevents tension cracks in plastic materials during sterilisation

### Application and dosage:

neodisher MediDry special is used in the final rinse in washer disinfectors. The dosing amount depends amongst other things on the water quality and the material of the items to be washed and is 0.3 - 1 ml/l.

Use suitable dosing devices.

### General notes on application:

- For professional use only.
- Do not mix with other products.

- Rinse out dosing system including suction hose with water before changing product.
- The reprocessing of medical devices should comply with all ordinances pursuant to the Medical Device Directive and should be performed with appropriate validated processes.
- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664.
- The instructions given by the manufacturer of the washer disinfectant are to be observed.

### Technical data:

pH-range	3.8 - 3.5 (0.3 - 1 ml/l, determined in deionised water, 20 °C)
Density	approx. 1.04 g/cm <sup>3</sup> (20 °C)
Viscosity	< 10 mPa s (concentrate, 20 °C)

### Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:

< 5% polycarboxylates

also preservatives

(methylchloroisothiazolinone/methylisothiazolinone)

### CE-mark:

neodisher MediDry special complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

### Storage information:

Always store at a temperature between 0 °C and 30 °C.

Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol



# neodisher<sup>®</sup> MediDry special

## Hazard and precautionary statements:

For further safety information, see EC safety data sheets. These are available e.g. at [www.drweigert.com](http://www.drweigert.com) under the category "Service/Downloads".

Dispose only when container is empty and closed.  
For disposal of product residues, refer to the Safety Data Sheet.

MB 4571/3-1 Revision date: 06/2020



With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.