





# neodisher® ProZyme

Product only available in the UK



## pH-neutral detergent for reprocessing thermostable and thermolabile instruments

Liquid concentrate

### Fields of application:

- Automated cleaning of thermostable and thermolabile instruments including surgical instruments, anaesthesia equipment, flexible endoscopes, containers and other medical equipment
- Manual cleaning of thermostable and thermolabile instruments in immersion baths or ultrasonic baths

### Performance spectrum:

- Removes residues of dried and denatured blood
- Suitable for stainless steel, instrument steel, light metals, glass, commonly used plastics and materials of anaesthesia equipment and flexible endoscopes

### Special properties:

- pH-neutral cleaning agent with enzymes and other cleaning adjuvants
- Excellent material compatibility

### Application and dosage:

neodisher ProZyme can be used in washer disinfectors as well as in immersion and ultrasonic baths. The dosing amount depends among other things on the respective field of application and the individual degree of contamination of the instruments to be reprocessed. The following parameters are recommended when using neodisher ProZyme:

Application recommendations	
Automated cleaning of thermostable and thermolabile instruments	2 - 10 ml/l (0.2 - 1.0%), 40 - 50 °C, 5 - 10 min
Automated cleaning of flexible endoscopes	5 ml/l (0.5%), 5 - 10 min
Manual cleaning of thermostable and thermolabile instruments in immersion baths	5 - 20 ml/l (0.5 - 2.0%), 20 - 30 min, max. 45°C

Manual cleaning of thermostable and thermolabile instruments in ultrasonic baths	5 - 20 ml/l, (0.5 – 2.0%) 5 - 10 min, max. 45°C
--	--

The working solution of neodisher ProZyme has to be replaced at least daily; immediately in case of visible soiling.

The neodisher ProZyme solution has to be rinsed off completely (preferably with deionised water). To avoid water stains, the use of fully demineralised water in the final rinse is recommended.

### Notes on application:

- For professional use only.
- Rinse out dosing system including suction hose with water before changing product
- For economical and controlled dosing we recommend using manual dosing ancillaries, such as e.g. mixing equipment or dosing caps. Please contact us.
- Reprocessing should comply with all ordinances pursuant to the Medical Device Directive and should be performed with appropriate validated processes.
- Please observe the medical device manufacturer's recommendations for use according to the requirements of DIN EN ISO 17664.
- Do not mix with other products.
- The instructions given by the manufacturer of the washer disinfectant are to be observed.

### Technical data:

pH value	7.9 (2 - 20 ml/l, determined in deionised water, 20 °C)
Viscosity	< 10 mPa s (concentrate, 20 °C)
Density	approx. 1.0 g/cm <sup>3</sup> (20 °C)



# neodisher<sup>®</sup> ProZyme

## Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:  
< 5 % anionic surfactants, also enzymes,  
preservatives (2-phenoxyethanol)


## CE-mark:

neodisher ProZyme complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

## Storage information:

Always store at a temperature between 0 °C and 25 °C. Keep away from sunlight.

Usable for 2 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

## Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available e.g. at [www.drweigert.com](http://www.drweigert.com) under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

MB 4756/3-4

Date of issue: 05/2021

---

With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.

---