



neodisher[®] SBN plus



Special rinse aid for use
in washer-disinfectors for bedpans



Liquid concentrate

Fields of application:

- Complexing water hardness and rinsing when reprocessing bedpans and urine bottles in bedpan washer-disinfectors with steam disinfection

Performance spectrum:

- Effectively prevents lime residues in steam generators, pipes and jets also at high degrees of water hardness
- With special components for effectively complexing water hardness
- With surfactants in a low-foam formula for optimized, rapid and stain-free drying without limescale deposits

Special properties:

- With a fresh and pleasant odour, reliably covers unpleasant odours
- Suitable for all degrees of water hardness
- Acidic, based on organic acids

Application and dosage:

neodisher SBN plus is used in bedpan washer-disinfectors with steam disinfection. neodisher SBN plus is dosed using dosing equipment installed in the washer-disinfector. The dosing amount is 0.5 - 1.5 ml/l depending on water hardness. Dosing systems and steam generators must be acid-compatible.

For removing stubborn soiling we recommend the additional use of neodisher SBR extra, a special detergent for bedpans and urine bottles.

Notes on application:

- For professional use only.
- Rinse out dosing system including suction hose with water before changing product.
- Reprocessing should comply with all ordinances pursuant to the regulations on

medical devices and should be performed with appropriate validated processes.

- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664.
- The instructions given by the manufacturer of the washer disinfector are to be observed.
- Do not mix with other products.

Technical data:

pH-range	3.5 - 3.1 (0.5 - 1.5 ml/l, determined in deionised water, 20 °C)
Viscosity	< 50 mPa s (concentrate, 20 °C)
Density	approx. 1.1 g/cm ³ (20 °C)

Ingredients:


5 - 15 % non-ionic surfactants
< 5 % polycarboxylates
also preservatives
(methylchloroisothiazolinone/methylisothiazolinone),
perfumes

CE-marks:

neodisher SBN plus complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

Storage information:

Always store at a temperature between 0 °C and 30 °C. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

Hazard and precautionary statements:

For safety information see Safety Data Sheets.



neodisher[®] SBN plus

These are available e.g. at www.drweigert.com
under the category "Service/Downloads".

Dispose only when container is empty and
closed. For disposal of product residues, refer to
the Safety Data Sheet.

MB 4792/3-2 Date of issue: 03/2020



With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.
