




neodisher[®] SBX



Complexing agent for the automated reprocessing of bedpans

Liquid concentrate

Fields of application:

- Complexing water hardness to prevent the formation of lime deposits during the automated reprocessing of bedpans and urine bottles

Performance spectrum:

- Prevents lime residues in steam generators, pipes and jets

Special properties:

- For sequestration of water hardness
- Suitable for all degrees of water hardness
- Mildly alkaline formulation

Application and dosage:

neodisher SBX is used in bedpan washer disinfectors with steam disinfection. neodisher SBX is dosed via the dosing system integrated in the bedpan washer disinfectant. The dosing amount is 1-3 ml/l depending on the degree of water hardness.

neodisher SBX can be dosed directly into the water used for steam generation or alternatively it may be dosed into a water reservoir in the bedpan washer disinfectant, from which the water is taken for the cleaning step and the steam generation.

For removing stubborn soiling we recommend the additional use of neodisher SBR extra, a special detergent for the automated cleaning of bedpans and urine bottles.

Notes on application:

- For professional use only.
- Rinse out dosing system including suction hose with water before changing product.
- Reprocessing should comply with all ordinances pursuant to the regulations on

medical devices and should be performed with appropriate validated processes.

- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664.
- Do not mix with other products.
- The instructions given by the manufacturer of the washer disinfectant are to be observed.

Technical data:

pH-range	10.1 - 10.2 (1 - 3 ml/l, determined in deionised water, 20 °C)
Viscosity	< 50 mPa s (concentrate, 20 °C)
Density	approx. 1.2 g/cm ³ (20 °C)

Ingredients:


Ingredients according to Regulation (EC) No 648/2004 on detergents:
15 - 30 % polycarboxylates

CE-marks:

neodisher SBX complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

Storage information:

Store in a cool and frost-free place. Always store at a temperature between 0 °C and 30 °C. Keep away from sunlight. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

Changes in the colour of the product may occur when storing in factory-sealed trade



neodisher[®] SBX

units. This has no impact on the properties of the product which are relevant for application.

Hazard and precautionary statements:

neodisher SBX is not a hazardous product according to the CLP-Regulation (EC) No 1272/2008.

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

For safety information see Safety Data Sheets. These are available e.g. at www.drweigert.com under the category "Service/Downloads".

MB 4015/3-4 Date of issue: 04/2024

With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.
