




neodisher[®] Septo DN



Disinfectant for the automated reprocessing of flexible endoscopes and thermolabile instruments

Liquid concentrate

Fields of application:

- High-level disinfection of flexible endoscopes and thermolabile instruments in washer disinfectors

Performance Spectrum:

- Bactericidal, fungicidal, mycobactericidal and virucidal activity. Confirmed by certification according to European Standards
- In the cleaning step the detergent neodisher MediClean forte is recommended
- The automated reprocessing with neodisher MediClean forte and neodisher Septo DN fulfils the requirements of DIN EN ISO 15883-4 regarding a germ reduction of > 9 log steps in the entire process.

Special Properties:

- Very low odour
- High material compatibility. Suitable for the endoscopes of all leading manufacturers
- Free of formaldehyde and quaternary ammonium compounds
- The disinfecting efficacy has been tested and confirmed according to EN 14885. Therefore, neodisher Septo DN complies with the requirements for disinfectants for disinfecting flexible endoscopes according to DIN EN ISO 15883-4
- On the IHO¹ list of disinfectants

Application and dosage:

neodisher Septo DN is used in washer disinfectors. neodisher Septo DN is dosed via the dosing device integrated in the washer disinfectant at the beginning of the disinfecting step. The following parameters must be observed to achieve the activity spectrum mentioned below:

bactericidal, fungicidal, mycobactericidal and virucidal activity	10 ml/l (1.0 %), 55 °C, 5 min
---	-------------------------------

For avoiding water stains the use of deionised water in the final rinse is recommended.

General instructions on use:

- For professional use only.
- The neodisher Septo DN working solution has to be rinsed off completely with water (preferably deionised water).
- Rinse out dosing system including suction hose with water before changing product.
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664.
- The instructions given by the manufacturer of the washer disinfectant are to be observed.
- Do not mix with other products.

Expert Reports:

The disinfecting activity has been confirmed by certification. Expert reports are available on request.

Technical data:

pH-value	4.3 (1.0 %, determined in deionised water, 20 °C)
Viscosity	< 10 mPa s (concentrate, 20 °C)
Density	approx. 1.0 g/cm ³ (20 °C)

Ingredients:

Active substances in 100 g:
10.5 g glutaral




neodisher Septo DN complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.



neodisher[®] Septo DN

Storage information:

Always store at a temperature between 0 °C and 30 °C. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available at www.drweigert.com under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

MB 4010/3-5
Revision Date: 06/2020

¹Industrieverband Hygiene und Oberflächenschutz [German Association for Hygiene and Surface Protection]



With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.
