



neiform MED rapid



Rapid disinfectant for medical devices

Ready-for-use solution



Fields of application:

- Specific and rapid disinfection of surfaces of medical devices such as surgical tables in hospitals as well as in doctors' and dentists' practices

Performance spectrum:

- Bactericidal, yeasticidal, fungicidal and tuberculocidal activity confirmed in accordance with EN methods¹ as well as activity against enveloped viruses confirmed according to the test guidelines of the RKI/DVV² by certification
- Also limited spectrum virucidal activity (active against noro and adenovirus) and active against rotavirus
- The surfaces to be disinfected must be compatible with alcoholic disinfectants
- Not suitable for disinfecting acrylic glass (perspex)
- Included in the IHO³ list of disinfectants

Special properties:

- Wide activity spectrum, even with short contact times
- Rapid and stain-free drying, no rinsing with fresh water necessary
- Free of aldehydes, colourants and perfumes

Application and dosage:

Apply neiform MED rapid undiluted on clean and dry surfaces and disinfect by wiping. Surfaces must be completely wetted with neiform MED rapid. To achieve the desired activity, the below-mentioned contact times must be strictly adhered to.

	Clean conditions, 20 °C (cleaned surfaces)
bactericidal (EN 13727, EN 16615)	undiluted, 1 min
yeasticidal (EN 13624, EN 16615)	undiluted, 1 min

fungicidal (EN 13624)	undiluted, 20 min
tuberculocidal (EN 14348)	undiluted, 30 sec
Active against enveloped viruses (incl. HIV, HBV, HCV) (EN 14476, RKI/DVV ²)	undiluted, 15 sec
Limited virucidal activity (adenovirus, norovirus, incl. all enveloped viruses) (EN 14476)	undiluted, 2 min
Active against norovirus (MNV ⁴) (EN 14476)	undiluted, 2 min
Active against adenovirus (EN 14476)	undiluted, 1 min
Active against rotavirus (EN 14476)	undiluted, 1 min

Due to its alcohol content, neiform MED rapid must not be applied in too large amounts (apply max. 50 ml per m² on the surfaces to be treated and max. 100 ml per m² base area). Please thoroughly read and observe the hazard and precautionary statements.

General instructions for use:

- For professional use only.
- For economical and controlled dosage the use of manual dosing aids is recommended. Please contact us.
- It is generally recommended to wear gloves when working with disinfectants.
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Please observe the reprocessing recommendations given by the manufacturer of the medical devices and surfaces. Do not mix with other products.



neoform MED rapid

Expert reports:

The disinfecting activity has been confirmed by certification.

Technical data:

Viscosity	< 10 mPa s (undiluted, 20 °C)
Density	approx. 0.9 g/cm ³ (20 °C)
Flash point	27 °C according to DIN 22791

Ingredients:


Active substances in 100 g:
35.0 g propan-1-ol,
25.0 g ethanol



neoform MED rapid complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

Storage information:

Always store at a temperature between -15 °C and 25 °C. Keep container tightly closed. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available at www.drweigert.com under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to Safety Data Sheet.

MB 4070/3-5 Date of Issue: 01/2024

- 1 European test methods
- 2 according to the guidelines of the Robert Koch-Institute (RKI) and the Deutschen Vereinigung zur Bekämpfung von Viruskrankheiten (DVV) [German Association for the Control of Virus Diseases]
- 3 Industrieverband Hygiene und Oberflächenschutz [Industry Association Surface Protection and Hygiene, Germany]
- 4 Murine norovirus

With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.