



# neoform<sup>®</sup> Classic



## Disinfecting cleaner for surfaces



### Liquid concentrate

#### Fields of application:

- Cleaning and disinfecting surfaces of medical devices, medical equipment and other surfaces in hospitals, doctors' practices, old people's homes, laboratories and public institutions
- Also suitable for use with the Dr. Weigert fleece wipe dispensing system neoform wipes RTF

#### Performance spectrum:

- Bactericidal, yeasticidal, tuberculocidal activity as well as activity against enveloped viruses and limited spectrum virucidal activity confirmed and certified
- Activity confirmed in the 4-field-test (EN 16615)
- Disinfecting activity has been tested and proved according to European standards and VAH<sup>1</sup> methods
- VAH<sup>1</sup>-listed, on the IHO<sup>2</sup> disinfectant list
- The application in the Dr. Weigert fleece wipe dispensing system neoform wipes RTF fulfils the requirements of the current VAH recommendation<sup>3</sup>

#### Special Properties:

- Wide activity spectrum, even at low concentrations resp. short contact times
- Good cleaning performance
- Free of aldehydes
- Excellent material compatibility

#### Application and dosage:

neoform Classic is used for the combined cleaning and disinfection of surfaces at room temperature via wiping. For achieving the desired activity dilute neoform Classic according to the below-mentioned dosing recommendations. Wipe the surfaces with the working solution. Adhere to the below-mentioned contact times.

For instructions on how to use the product with the neoform wipes RTF fleece wipe

dispensing system please refer to the neoform wipes RTF data sheet. The soaked neoform wipes can be used for up to 28 days when the wipe system is kept tightly closed.

Disinfecting performance (20°C) dirty conditions (non-pre-cleaned surfaces)	
Bactericidal VAH, EN 13727, EN 16615)	5 ml/l (0.5 %), 5 min
yeasticidal (VAH, EN 13624, EN 16615)	5 ml/l (0.5 %), 5 min
tuberculocidal (EN 14348)	15 ml/l (1.5 %), 60 min
activity against enveloped viruses (incl. HIV, HBV, HCV) (EN 14476)	10 ml/l (1.0 %), 15 min 5 ml/l (0.5 %), 30 min
limited spectrum virucidal activity (norovirus, adenovirus, incl. all enveloped viruses) (EN 14476)	20 ml/l (2.0 %), 120 min
Rotavirus (EN 14476)	20 ml/l (2.0 %), 90 min

When used on sensitive surfaces of medical devices, we recommend wiping with water of at least drinking water quality after the contact time has elapsed.

Aldehydic products should not be used before or after disinfection with neoform Classic. If an aldehyde containing product was used beforehand an intermediate rinse should be carried out.

#### General notes on application:

- For professional use only!
- For the economical and controlled dosage, the use of manual dosing aids such as mixing equipment or dosing caps is recommended. Please contact us.



# neiform<sup>®</sup> Classic

- It is generally recommended to wear gloves when working with disinfectants.
- Processing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Do not mix with other products.
- *Use disinfectants safely. Always read the label and product information before use.*

## Expert reports:

The disinfecting activity has been tested and confirmed by certification. We will be pleased to provide certificates on request.

## Technical data:

pH-value	10 (concentrate, 20 °C)
Viscosity	< 10 mPa s
Density	approx. 1.0 g/cm <sup>3</sup> (20 °C)

## Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:

Disinfectant

active substances in 100 g:

9 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine

0.9 g amines, N-C12-14-(even-numbered)-alkyltrimethylene derivates, reaction products with chloroacetic acid




neiform Classic complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

## Storage information:

Always store at a temperature between 0 °C and 30 °C. Usable for 3 years when stored as recommended. For expiry date refer

to the stamp mark on the label behind the hourglass symbol .

Changes in the colour of the product may occur when storing in factory-sealed trade units. This has no impact on the properties of the product which are relevant for application.

## Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available at [www.drweigert.com](http://www.drweigert.com) under the category "Service/Downloads".  
Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

MB 4470/3-6

Revision date 07/2024

<sup>1</sup> Verbund für Angewandte Hygiene [Association of Applied Hygiene, Germany]  
<sup>2</sup> Industrieverband Hygiene und Oberflächenschutz [Industry Association Surface Protection and Hygiene, Germany]  
<sup>3</sup> Empfehlung zur Kontrolle kritischer Punkte bei der Anwendung von Tuchspendensystemen im Vortränksystem für die Flächendesinfektion (Hyg Med 2012; 37-11) [Recommendation for checking critical points when using wipes dispenser systems in a pre-soak system for surface disinfection]

With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.